

PARticle Therapy REsearch Center (PARTREC) From Nuclear Physics to Medicine

Prof. Alexander Gerbershagen

Head of Accelerator Facility

a.gerbershagen@umcg.nl







- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC
 - Nuclear Physics
 - Nuclear and Molecular Medicine
 - Radiation Biology
 - Clinical Translation Examples
 - Minibeams
 - Microdosimetry
 - Fast Irradiations
 - FLASH
 - VHEE
 - Novel Gantries
 - Patient Imaging




University Medical Center Groningen

















- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC
 - Nuclear Physics
 - Nuclear and Molecular Medicine
 - Radiation Biology
 - Clinical Translation Examples
 - Minibeams
 - Microdosimetry
 - Fast Irradiations
 - FLASH
 - VHEE
 - Novel Gantries
 - Patient Imaging




University Medical Center Groningen

partrec












PartrecGroningen University

- Rijksuniversiteit Groningen (RUG):
 - Founded in 1614 (second oldest in the Netherlands)
 - Ranked 76th in the world, 2nd in the Netherlands (Shanghai University Ranking)
 - Discoveries/Inventions
 - Calculus (Bernoulli 1720's and later)
 - Electric cars (Stratingh, 1837)
 - Vitamine C (Szent-Györgyi, 1928)
 - Phase contrast microscopy (Zernike, 1930s) and much more!
 - 4 Nobel Prize winners, incl. Feringa (Chemistry, 2016)
- KVI/PARTREC
 - Only accelerator physics institute in NL
 - Reorganization, inclusion of medical and radiobiological research
 - -> Re-established as PARticle Therapy REsearch Center (PARTREC)

















University Medical Center Groningen

Integrated with Medical Faculty of University of Groningen



partrec

Integrated into UMCG in 2019








First patient in 2018






PartrecAGOR Cyclotron

- Superconducting AGOR cyclotron is a multi-particle, variable energy AVF-cyclotron
- French-Dutch collaboration built 1987 1994
- Operational since 1996
- Magnetic field (1.7 to 4.1 T) produced by
 - Two pairs of superconducting main coils
 - fifteen trim coils
 - three iron hill sectors for focussing
- 3 halfwave RF cavities, 24 62 MHz; h = 2, 3 or 4
- Three external ion sources (two ECR sources for protons and heavy ions, multi-cusp source for light ions) are axially injected
- Extraction
 - 300 500 turns depending on harmonic mode
 - extraction radius 870 890 mm depending on E/A
 - turn separation at extraction 2 3 mm ~ beamwidth






Beams and Test Rooms



AGOR can deliver beams of all elements up to Xe









partrecOur Team

- Faculty (4)
- Post-doc, PhD-students (~10)
- Technical staff (25)
 - Operators
 - Experimental and project support
 - Cryogenics, cooling and vacuum
 - Design/Mechanical
 - Electronics
 - IT Support
- Operational 120 hours/week, 26 weeks/year
- Users: Medical Physics, Radiobiology, ESA, Detectors, Commercial
- Beam requests: <u>irradiations.partrec@umcg.nl</u>








IMPACT Beam Line

- >3 MEUR grant from KWF
- Designed for small animals
- New beam line
 - Pencil beam scanning (10 cm x 10 cm)
 - Shoot-through or SOBP
 - Collimator for penumbra
- Combined on-line 3D X-ray imager and X-ray irradiator
- Irradiation planning
- First beam: 04.10.2024
- Fully operational: 2025










PartrecOne Stop Shop

- Experiment development
- Ethics authorisation process
- Animal procurement logistics
- On site animal accommodation with IVCs
 - capacity 200 rats and mice
 - no long term stay
 - two additional accommodations planned
- "Twin Beam" Irradiation
- Laboratory for animal handling prior and post irradiation
 - GronSAI imaging center: optical, molecular, CT, MRI, PET
- Data management facilities









- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC





- Nuclear Physics
- Nuclear and Molecular Medicine
- Radiation Biology
 - Clinical Translation Examples
 - Minibeams
- Microdosimetry
- Fast Irradiations
 - FLASH
 - VHEE
- Novel Gantries
- Patient Imaging




University Medical Center Groningen













PartrecHeavy Ion Beams


- AGOR can deliver beams of all elements up to Pb
- Research areas
 - Radiobiology (RuG, UMCG, PSI)
 - Detector tests & development (ESA)
 - Experiment development (ESA)
 - Radiation hardness (ESA, companies)

NEXT (J. Even)

- New experimental research on the production of neutron-rich heavy nuclei using multi-nucleon transfer reactions between heavy nuclei (e.g. ¹³⁶Xe on ²⁰⁸Pb) has recently been started
- ECR ion source development, improvement of transmission from source to extraction
- A new experimental station consisting of a 3 T superconducting solenoid fragment separator and MR-ToF mass spectrometer is developed with RUG has been installed









- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC
 - Nuclear Physics
- Nuclear and Molecular Medicine
- Radiation Biology
 - Clinical Translation Examples
 - Minibeams
- Microdosimetry
- Fast Irradiations
 - FLASH
 - VHEE
- Novel Gantries
- Patient Imaging




University Medical Center Groningen

partrec













partrec **Terbium Theragnostics**

- >10 MEUR grant from IPCEI
- Produce four Terbium isotopes for patient diagnostics and treatment
- Partners: SHINE, UMCG



Tb 149

4.1 h


y 352:

Tb 155

5.32 d

180, 262

4.2 m




The concept of theragnostic radionuclide pairs


(a) Diagnosis



partrec NuCapCure

- Consortium of eight parties
- Goal: Treat Glioblastoma (presently <10% 5-year survival)
- Envisaged radiation effect:
 - Higher in tumour due to boron proton capture
 - Lower in healthy tissue due to FLASH
- > 5 MEUR from EIC Pathfinder Open 2023
- Ranked 1st (788 proposals submitted, 56 funded)










University Medical Center Groningen


Image-Guided Pharmacotherapy

IMAGING REVEALS LOCATION



RADIATION ACTIVATES THERAPY





Radiation setups, molecules, activation modes

- ➤ 18 MEUR grant from HTRIC x UEF
- Proposal headed by Prof. Ben Feringa, Nobel prize laureate







- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC
 - Nuclear Physics
 - Nuclear and Molecular Medicine
 - Radiation Biology
- Clinical Translation Examples
- Minibeams
- Microdosimetry
- Fast Irradiations
 - FLASH
 - VHEE
- Novel Gantries
- Patient Imaging




University Medical Center Groningen

partrec















Stem cell sparing radiotherapy

Rat model



Patients



Stem cell region

- Centre of the glands is rich in organ-specific stem cells.
- Irradiation of this sub-volume results in degeneration of the entire gland.
- Irradiation of the other parts only causes local damage
- Retrospective analysis of a Canadian patient cohort
 - dose to the centre of the gland is the best predictor of post-treatment saliva production.
- Validated in a trial at UMCG.
- Currently implemented as the standard approach to minimize the risk of radiation-induced xerostomia in patients






Steenbakkers et al. IJROB 2022



Heart/lung → Cardiopulmonary





- Van Luijk et al. Canc. Res. 2005
- Ghobadi et al. IJROB 2012
- Ghobadi et al. Thorax 2012



- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC
 - Nuclear Physics
 - Nuclear and Molecular Medicine
 - Radiation Biology
 - Clinical Translation Examples





- Minibeams
- Microdosimetry
- Fast Irradiations
 - FLASH
 - VHEE
- Novel Gantries
- Patient Imaging




University Medical Center Groningen

partrec












There is more than volume...




Effects of:

- Shape of dose distribution
- Regional differences
- Interaction between organs

Inherent loss of spatial data:

DVH does not tell location of a high or low dose region:

DVH does not tell how certain dose regions are distributed:



A single hot spot


Many small hot spots






PartrecSerial and Parallel Organs


- What and how to spare?
 - Many degrees of freedom:
 - Reduce high-dose regions? (many beams)
 - Reduce low-dose volume? (fewer beams)
 - Both? (Need costly proton/ion therapy)
- Responses are tissue/organ dependent!
- ⇒ Serial vs parallel organs
- Practical approximation used clinically:
 - Serial → Maximum dose as predictor
 - Parallel → Mean dose as predictor



PartrecSpinal cord

Homogeneous irradiation of different lengths of spinal cord












Spinal cord: inhomogeneous dose distribution







Spinal cord: inhomogeneous dose distribution




(1x 4 mm: ED50=53 Gy)

Volume-Effect? Configuration determines outcome






Spinal cord: include repair in model










Spinal cord: include repair in model

Correct the dose distribution for repair from unirradiated tissue





van Luijk et al. Int J Radiat Oncol Biol Phys. 2005

- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC
 - Nuclear Physics
 - Nuclear and Molecular Medicine
 - Radiation Biology
 - Clinical Translation Examples
 - Minibeams




- Microdosimetry
- Fast Irradiations
 - FLASH
 - VHEE
- Novel Gantries
- Patient Imaging




University Medical Center Groningen
















Tumour survival and cell killing




Irradiation with X-Rays

- What went wrong here?
- What's the relation between cell death and probability of curing a tumour?

- Low probability of dying from specific hit
- High number of hits
- → Poisson statistics!
 - SF = e^{-m}
 - SF: Surviving fraction

early cell death (apoptosis, senescence autophagy, necrosis) em: mean hits per cell (prop. to dose)








partrecMicrodosimetry vs LET

- The Microdosimetric Kinetic model (MKM)
 uses microdosimetry to describe the radiation
 field within each voxel and allows the RBE to be
 calculated.
- Useful tool to assess the energy deposited taking into account the stochastic fluctuations occurring at the nucleus scale
- Provides a more advanced description of radiation quality because:
 - Measures the energy deposited in a volume comparable to the cells nucleus
 - It takes into account stochastic fluctuations in energy deposition



Geant4:

- + widely validated for microdosimetry and clinical application
- + can score individual electrons' energy deposit with an arbitrary accuracy
- modelling a specific patient is very difficult
- illegal to use for treatment planning

Raystation:

- + easy patient modelling from CT scan
- + very fast proton, neutron, and alpha transport
- proton energy deposit modelled via tabulated stopping power
- spatial resolution ≥ mm





- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC
 - Nuclear Physics
 - Nuclear and Molecular Medicine
 - Radiation Biology
 - Clinical Translation Examples
 - Minibeams
 - Microdosimetry





- Fast Irradiations
 - FLASH
 - VHEE
- Novel Gantries
- Patient Imaging




University Medical Center Groningen















Mitigating organ / tumor motion

Organ motion




Possible solutions:

Gating



 Adaptive scanning (tumor tracking)

Fast rescanning





PartrecFLASH Therapy

- FLASH is a biological phenomenon and not defined via specific physical beam parameters
- Healthy tissue effect (not tumour tissue!)
- What have been the experimental conditions to observe a FLASH effect?
 - Small volumes of normal tissues (a few cc)
 - Mainly with single dose (> 7-8 Gy)
 - Generally with Overall Treatment Time (OTT) < 200 ms

25 Gy no necrosis

28 Gy

31 Gy

34 Gy

Conventional

FLASH



28 Gy

31 Gy


34 Gy, no necrosis






FLASH Irradiations

- Can deliver proton and helium beams with FLASH intensity
- Establish Twin Beams, replicating dose delivery parameters and control methods in conformity with clinical facilities




Groningen Proton Therapy Center (GPTC)


university of
groningen



Beam Dosimetry at FLASH Dose Rates

- FLASH therapy requires >40 Gy/s
- For human tumours (~1 liter) that means 100s nA of current
- Standard solution: ionization chambers saturate
- Proposed solution: Scintillating gas and camera





2D dose-monitor



PartrecVery High Energy Electrons

 Conventional proton or ion irradiations are performed with Spread-Out-Bragg-Peak (SOBP)


No technical solution to generating SOBP with same


dose conformity at FLASH time scale


=> Applied in "transmission mode"

Possible solution: electrons with > 100 MeV

- VHEE are less expensive
- VHEE are more sustainable
 - lower energy consumption
 - less material needed (iron, steel, concrete)
 - less activated materials














Dose Delivery for VHEE


- Scanning at FLASH time scale (100 ms) is impractical
- Solution: Scattering
- Tests at CLEAR successful
 - => Now a standard tool for dose distribution for radiobiological experiments at CLEAR





Thin initial beam





- In collaboration with
 - M. Dosanjh, C. Robertson (Oxford)
 - R. Corsini, A. Latina (CERN)









Beam magnified by flat first

scatterer (S1)

Beam flattened by Gaussian

second scatterer (S2)

Beam Monitor for VHEE


- Two arrays of silicon fibre sensors
- Measure Cherenkov radiation (not scintillation)
- Readout: CMOS camera
- Measurements at CLEAR
 - Profile measurements
 - Linearity of response with dose rate (up to 40 Gy per 100 ns pulse)
- In collaboration with


university of groningen


- M. Dosanjh, J. Bateman (Oxford)
- R. Corsini, I. Ortega (CERN)
















- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC
 - Nuclear Physics
 - Nuclear and Molecular Medicine
 - Radiation Biology
 - Clinical Translation Examples
 - Minibeams
 - Microdosimetry
 - Fast Irradiations
 - FLASH
 - VHEE
 - Novel Gantries
 - Patient Imaging







University Medical Center Groningen














Partrec Conical gantry - Commercial standard layout



IBA Sumitomo Hitachi Mitsubishi Varian

First commercial scanning gantry of Varian in Munich






- 135° bending magnet
 - Shorter length but larger radius
 - · Cylindrical treatment cell
- · Initially only for passive scattering
- Lately also for scanning






partrecGaToroid

- Static gantry
- Irradiation angle set by vector magnet
- Beam focussed and deflected via toroidal magnets
- Advantages: extremely fast stereotactic irradiation
 => Optimal for FLASH, motion mitigation, fast-Arc therapy
- Four versions:
 - Proton, ion, VHEE and proton pre-clinical
- Applying for funding to install at PARTREC
- In collaboration with University of Oxford and CERN
 (L. Bottura, M. Dosanjh, A. Haziot, A. Latina, T. Lehtinen, C. Robertson)











partrec Arc gantry

- **Outer Coils:** Generate magnetic field
- Inner Coils: Create a magnetic free area for treatment
- ISM (Injection Scanning Magnets): Set the treatment angle
- Challenges: Beam focussing (independent of treatment angle), scanning
- In collaboration with:















- Introduction:
 - University of Groningen, UMCG and PARTREC
 - Our accelerator
 - Our team, services and users
- Research at PARTREC
 - Nuclear Physics
 - Nuclear and Molecular Medicine
 - Radiation Biology
 - Clinical Translation Examples
 - Minibeams
 - Microdosimetry
 - Fast Irradiations
 - FLASH
 - VHEE
 - Novel Gantries
- Patient Imaging







University Medical Center Groningen













Real-time In-vivo VERification of proton therapy (RIVER)

Aim: "reducing treatment planning safety margins without reducing treatment safety"

How?

- PET imaging of the short-lived isotope N-12 ($T_{1/2} = 11 \text{ ms}$)
- Range determination per pencil beam spot


Patient benefit


- reduced safety margins lead to
- less dose to organs-at-risk
 - less complications
 - better quality of life

Research focus


- experiments with a head phantom
- mimic clinical workflow from TPS via Monte Carlo prediction to comparison with experimental PET images


experiments at PARTREC





calculations/simulations










Proton Radiography

- Allows precise tissue density measurement
- Two different methods:
 - Integral charge yield (at single or multiple energies)
 - Individual particle tracking
- Irradiations of phantoms and mice performed at **PARTREC**


(LANEX screen + CCD Camera, 50 mm proton field)

Worldwide first clinical implementation at GPTC!



CBCT image generation

- fastCAT CBCT simulator
- MOBY/ROBY phantoms




CT conversion




Treatment planning

Dose calculation

- beam modelling w/ BDSIM
- irradiation plan from matRad
- MC dose calculation in TOPAS









- SECT, DECT
- Validation of calibration by proton imaging



partrec Summary

- UMCG has unique combination of treatment facility (GPTC) and research center accelerator facility (PARTREC)
- PARTREC delivers protons and ions (He, C, O and others)
- PARTREC has team of 40 people
- Ongoing research topics are:

Acknowledgement

Research funding

Access funding

 Radiation Biology, Microdosimetry, Minibeams, Ultra-fast Irradiations, Gantries, Patient Imaging, Nuclear Medicine

Host institutions and all colleagues for contributing

DUTCH CANCER



Netherlands Organisation for Scientific Research





University Medical Center Groningen

















